

# Tubular Reactor Systems



Series Number:

## Series 5400 Continuous Flow Tubular Reactor Systems

**Bench Top or** Floor Stand

Vessel Sizes, mL: 5 mL - 1000 mL

Standard Pressure MAWP Rating, psi (bar):

1500 (103) 3000 (207) 4500 (310)

Maximum Operating Temperature, °C:

350 or 550



Model 5403 with a 1" inside diameter x 24" length, 3-zone split tube furnace with gas feed system, cooling condenser, and gas/liquid separator vessel.

Tubular reactors are used in a continuous flow mode with reagents flowing in and products being removed. They can be the simplest of all reactor designs. Tubular reactors are often referred to by a variety of names:

- Pipe reactors
- Packed-bed reactors
- Fixed-bed reactors
- Trickle-bed reactors
- Bubble-column reactors
- Ebullating-bed reactors


Single-phase flow in a tubular reactor can be upward or downward. Two-phase flow can be co-current up-flow, counter-current (liquid down, gas up) or, most commonly, co-current down-flow.

Tubular reactors can have a single wall and be heated with an external electric furnace or they can be jacketed for heating or cooling with a circulating heat transfer fluid. External furnaces are typically rigid, split-tube heaters. Tubular reactors are used in a variety of industries:


- Petroleum
- Petrochemical
  - Polymer
- Pharmaceutical
- Waste Treatment
- Specialty Chemical
- Alternative Energy

Tubular reactors are used in a variety of applications, including:

- Carbonvlation
- Dehydrogenation
- Hydrogenation
- Hydrocracking
- Hydroformulation
- Oxidative decomposition
- Partial oxidation
- Polymerization
- Reforming



It is often desirable to size a tubular reactor to be large enough to fit 8 to 10 catalyst particles across the diameter and be at least 40-50 particle diameters long. The length to diameter ratio can be varied to study the effect of catalyst bed length by equipping the reactor with "spools" placed into the bottom of the reactor to change this ratio.





Open 3-zone Split Tube Furnace with 1" I.D. Tubular Reactor.

Tubular reactor systems are highly customizable and can be made to various lengths and diameters and engineered for various pressures and temperatures, and materials of construction.

A split-tube furnace is provided for heating these vessels. Insulation is provided at each end to minimize heat losses and prevent the end caps from being heated. The heater length is normally divided into one, two, or three separate heating zones, although it can be split into more zones if required.

A fixed internal thermocouple in each zone can be furnished or a single moveable thermocouple in a centerline thermowell can be used to measure the temperature at points along the catalyst bed. External thermocouples are typically provided for control of each zone of the heater, as can be seen in the photo above.

#### **Gas Feed Systems**

Various gas feeds can be set up and operated from a **Gas Distribution Rack**. In order to deliver a steady flow of gas to a reactor, it is necessary to provide gas at a constant pressure to an electronic **Mass Flow Controller**. This instrument will compare the actual flow rate delivered to the set point chosen by the user, and automatically adjust an integral control valve to assure a constant flow.

#### **Series 5400 Tubular Reactor System Specifications**

Shaded bar indicates specifications that change within series.

Model Number	5401	5402	5403	5404	
Sizes	3/8 in.	1/2 in.	1.0 in.	1.5 in.	
O.D. / I.D. (in.)	0.38 / 0.28	0.50 / 0.37	1.50 / 0.99	1.88 / 1.44	
O.D. / I.D. (mm)	9.5 / 7.0	13 / 9.5	38 / 25	48 / 36	
Heated Length (in.)	6, 12, 24		12, 24, 36		
Maximum Pressure (psi/bar)	3000 (207)		5000 (345)	3000 (207)	
Maximum Temperature	550 550 350		350		
Support Spools	N	0	Optional		
Spiral Pre-Heat	N	No Optional			
No. Ports in Top Head	1			4	
No. Ports in Bottom Head	1		4		
Internal Thermocouple	Optional (Moveable or multi-point fixed)				

Care must be taken to size these controllers for the specific gas, flow rate range, and maximum pressure of operation. A mass flow controller needs a power supply and read-out device, as well as a means of introducing the desired set point.

When ordering mass flow controllers, you will need to specify:


- 1. Type of gas to be metered (e.g. N<sub>2</sub>, H<sub>2</sub>, CH<sub>4</sub>)
- 2. Maximum operating pressure of the gas (100 or 300 bar)
- Maximum flow rate range in standard cc's per minute (sccm)
- 4. Pressure for calibration of the instrument


Mass flow controllers are available for use to 1500 psi and to 4500 psi. Considerable savings can be obtained if the mass flow controller is to be used only to 1500 psi.

The schematic at right depicts the installation of a mass flow controller for the introduction of gas to a continuous-flow reaction system. Such installations are enhanced with the addition of a by-pass valve for rapid filling or flushing.

A purge line can also be added. It is typically used for feeding nitrogen or helium to remove air before reaction or to remove reactive gases before opening the reactor at the end of a run. The purge line includes a shut-off valve, filter, metering valve, and a reverse-flow check valve.

Shut-off valves can be automated when using a 4871 Control system.





## Series 5400 Continuous Flow Tubular Reactor Systems



Up-flow 5403 Tubular Reactor System with 300 mL heated volume, one purge line, one gas feed, two liquid feeds, product cooling condenser, and automated 2-phase back pressure regulator. An automated liquid sampler captures representative samples at a user-programmable interval.

#### **Liquid Metering Pumps**

High pressure piston pumps are most often used to inject liquids into a pressurized reactor operating in a continuous-flow mode. For low flow rates, HPLC pumps, many of which are rated for 5000 psig, are excellent choices. Typical flow rates for pumps of this type range up to 10 or 40 mL per minute. Pumps are available to accommodate manual control from their digital faceplate or computer-control from a 4871 Process Controller.

Chemical feed pumps are our recommendation for continuous feeding of liquids when the desired flow rate is greater than 2 liters per hour. Parr can assist with the feed pump selection. We will need to know the type of liquid; the minimum, typical, and maximum desired feed rate; the maximum operating pressure; and any special operating considerations such as corrosion possibilities.

#### **Cooling Condensers**

Cooling condensers are available to cool the products of the reaction. An adaptation of our standard condensers provides an excellent

#### **Back Pressure Regulators**

The reactor pressure is maintained by a Back Pressure Regulator (BPR) installed downstream of the reactor. This style of regulator will release products only when the reactor pressure exceeds a value preset by the operator.

When a BPR is used in conjunction with mass flow controllers, the user can maintain a constant flow of gas through a reactor held at an elevated constant pressure. This provides for the highest degree of control and reproducibility in a continuous-flow reactor system.

An alternate BPR may be available to permit pressure let-down of a two-phase stream from the reactor. This BPR requires that the operator provide a source of nitrogen or air at a pressure slightly above the desired operating pressure. With this style of BPR, the highpressure gas/liquid separator can be replaced with a low-pressure liquid product receiver allowing collection of near real-time liquid samples.



This continuous flow stirred reactor system is on a cart with our Modular Frame System. This modular frame allows for easy access and flexibility in hook-ups, accessories, and flow, including an interchangeable tubular reactor.

#### **Gas/Liquid Separators**

Tubular reactors operating in continuous-flow mode with both gas and liquid products will typically require a Gas/Liquid Separator. The separator is placed downstream of the reactor, often separated from the reactor by a cooling condenser. In the separator vessel, liquids are condensed and collected in the bottom of the vessel. Gases and non-condensed vapors are allowed to leave the top of the vessel and pass to the back pressure regulator. It is important to operate the standard BPR with a single fluid phase to prevent oscillation of the reactor pressure.

The gas/liquid separator can be sized large enough to act as a liquid product receiver that is drained periodically. Many of the non-stirred pressure vessels made by Parr are ideally suited for use as gas/liquid separators. Vessels of 300, 600, 1000, or 2000 mL are commonly chosen. Upon request, the bottom of the separator can be tapered to facilitate draining.

### **Control and Data Acquisition Systems**

A variety of solutions exist to meet the needs of system operators. System accessories such as heaters, mass flow controllers, and pumps can be obtained with individual control packages to create a Local Control System (LCS) based on our 4838 and 4848 Controllers.

As the number of channels to be controlled increases, economics and convenience will often dictate that the system of individual controllers should be replaced with the computer-based Model 4871 Process Controller (PCC).




Down-flow 5402 Tubular Reactor System with 20 mL volume, automated shutoff valves to safely shutdown in case of alarm/interlock, two gas feeds, one purge line, one liquid feed, and automated back pressure regulator. Automatically alternates between high pressure liquid collection and draining, with one of two heated gas/liquid separators collecting while the other drains.


## Series 5400 Continuous Flow Tubular Reactor Systems

On this page are schematic representations of typical tubular reactor systems, along with a symbols chart to facilitate understanding. We have provided an ordering number for each of these examples.

#### **Key to Symbols**




Order No. for this system would be: **5402B-SS-115-ST1(6)-1500-LCS-GF(1)-PL-LF(1)-ITW-CCD-GLS(300)-MPC** 



Order No. for this system would be:

5403F-SS-230-ST3(24)-3000-PCC-GF(2)-PL-LF(1)-ISP-CSS-ITW-GLS(600)-APC-ASV(3)



## Series 5400 Ordering Guide

A composite identification number to be used when ordering a 5400 Series Reactor can be developed by combining individual symbols from the separate sections below.

A Base Model		
Model No.	Size (0.D. / I.D.)	
5401	3/8 in. (0.38" / 0.28")	
5402	1/2 in. (0.50" / 0.37")	
5403	1.0 in. (1.9" / 1.0")	
5404	1.5 in. (2.0" / 1.5")	

Add suffix F for Floor Stand mounting Add suffix B for Bench Top mounting

B Materials of Construction		
-SS	T316 Stainless Steel	
-HC	Alloy 276	
-TI	Titanium	
-IN	Alloy 600	
-M0	Alloy 400	

C Electrica	al Supply
-115	115 VAC, 50/60Hz
-230	230 VAC, 50/60Hz

D Heater C	Options
-ST1(#)	Split Tube, 1-zone
-ST3(#)	Split Tube, 3-zone
-WJ(#)	Welded Jacket

Add suffix (6), (12), (24), (36) for heated length (in.)

E Maximum Operating Pressure		
-1500	1500 psi / 103 bar	
-3000	3000 psi / 207 bar, 200 bar for CE	
-4500	4500 psi / 310 bar	

<b>(F)</b> Controller			
-PCC	PC-based Process Control (4871-style)		
-LCS	Local Control System (4838-style)		

<b>G</b> Custom	Options		
-GF(#)	Number of Gas Feeds		
-PL	Purge Gas Feed Line		
-LF(#)	Number of Liquid Feeds		
-ISP	Internal Pre-heat Spiral (5403/5404 only)		
-CSS	Catalyst Support Spools (5403/5404 only)		
-ITW	Internal Thermowell, with Moveable T/C		
-IZT	Internal, 3-PT, Fixed T/C		
-CCD	Cooling Condenser		
-GLS(#)	Gas/Liquid Separator (300, 600, 1000, 2000 mL)		
-SPH	Separator Heater		
-MPC	Manual Pressure Control		
-APC*	Automated Pressure Control		
-ASV(#)*	Automated Shut-off Valves (1-12)		
-ASV(#)*	Automated Shut-off Valves (1-12)		

*Available	only with	/1071	Droocc	Control	/DCC\
"Avallable	OHIV WHILE	40 / I	Process	Common	IPLLL

H Certifications		
-No Symbol	No Certification Required	
-PARR	Parr Certification	
-ASME	ASME Certification	
-PED	PED Certification	



The system above has three 250 mL tubular reactors operating in parallel and controlled by a 4871 Process Controller with operator interface on a single PC. This system has weighed and/or heated feed tanks and a two-stage pressure let down.