

# ECS8024 | ORGANIC ELEMENTAL ANALYZER NC SOIL SPECIAL


## ECS 8024 - NC SOIL SPECIAL Organic Elemental Analysis

ECS 8024 is a specialized N and C Elemental Analyzer Model based on the Dumas combustion method.

ECS 8024 is prepared and dedicated to the carbon and nitrogen determination in soils and similar samples. ECS 8024 is based on sample combustion and separation of gases with a chromatographic column.

The combustion products, i.e. CO2 and N2 are separated and quantified by a high resolution TCD detector.





ECS 8024 is fully automatic, from the available samplers to the oxygen dosage, from the monitoring of consumables status to the automatic leak test



## ECS 8024 is composed by:

Dual furnace combustion system (for a better combustion and optimization of catalysts)
Safety quick fit system for reactors connection (easy and safe way to set the instrument)
Customized GC column
Detection system
Data acquisition and handling

Several configurations may be set for the determination of the target element. Instrument is already prepared for soil analysis and carbon and nitrogen determination. ECS 8024 could be also equipped with other tools for the upgrade of the instrument to the CHNS-O analyzer.



Consumables kit are prepared for obtaining the best analysis performances.

Prepacked reactors and all chemicals and tools allow fast and user-friendly analyses.

The versatility of ECS 8024 can be expressed also as:

Sample weight Sample type Liquid or solid sample





## ECS 80 SERIES

## Three different samplers are available:

Pneumatic autosampler up to 147 positions Electronic autosampler with 32, 50, 100 positions Manual sampler

Automatisms make it particularly user-friendly:

Automatic oxygen dosage system

(for a better consumption of oxygen and consumables)

Automatic consumables status monitoring

Automatic leak test

Standby mode

(gas, energy and time savings)





Automatic oxygen dosage: main highlights.

Details of sample are registered:

Position

Type

Weight

ECS 8024 will dose the oxygen for the right combustion. No oxygen or consumables wastes



## **EAS Clarity**

Powerful software for powerful instrument Getting data, analyzing data, presenting data

#### Integration

There is extensive possibility to modify chromatograms. The chromatogram can be changed by entering global parameters or interactively, through direct graphic modification of the baseline.

#### Overlay


Simultaneously displays a virtually unlimited number of chromatograms and their mathematical modification.


#### Calibration

Internal and external standard calculation methods.

#### Automated measuring support

Sequence tables for any set of samples.







### ECS 80 SERIES

#### What the user can do

#### Summary result tables

Displays and prints selected results from all simultaneously displayed chromatograms.

#### User settings

User selects parameters for peak display and the specification for axes, including color from an extensive array of color settings. Text labels and lines, either as part of the area or anchored to a chromatogram, may also be inserted.

#### **Export**

Optional exportation of all results, with or without the chromatogram, in various formats, into a file or clipboard.

#### **Import**

Imports chromatograms or mathematical curves, which have been saved in text or AIA formats, from other programs.

#### Special features

#### LIMS

Clarity offers connectivity with LIMS both for sample submission and result output. This can be done via convenient ASCII transfers.

#### Method and calibration history

Each chromatogram can easily be displayed under the same conditions as when it was printed, exported or saved.

#### Column performance

Calculations of peaks in terms of symmetry, efficiency, resolution; all by several methods.

#### Batch

Automatically batch processes, displays, exports or prints any number of chromatograms.

#### Language localization available

Basic version in English language.

At the moment, there are French, German, Russian, Spanish and Chinese localizations available in single installation package















## ECS 80 SERIES

#### ECS 8024 key points

- √ Fully automated analysis system
- √ High sensitivity, accuracy and precision
- ✓ Flexibility and versatility of applications
- ✓ Detector does not require reference gas
- ✓ Powerful software for viewing results from a computer
- ✓ Touch-screen display for an easy settings management

- √ Consumables status monitoring for an optimization of catalysts usage
- ✓ Three types of samplers available: electronic / pneumatic / manual.
- ✓ Easy connection to Mass Spectrometers and other detectors for stable isotopes analysis
- √ Low operation and management costs
- ✓ Standby mode for gas, energy and time savings





#### ECS 8024 application fields

- ✓ Soil science and geology- marine science (distinction between organic and inorganic carbon is available through the previous acidification of the sample)
- ✓ Environmental analysis









## **Analytical and Technical Features**

## **ECS 8024 Features**

Type CN

Upgradable to CHNS-O

Analysis time CN 5 min

Analytical range C 0.002-20 mg

0.002-20 mg

Accuracy\* <0,2% (certified standard; purity >99.9%)
Precision\* <0,1% (certified standard; purity >99.9%)

, 10000.00.

Sampler Pneumatic autosampler 147 positions max Electronic autosampler 32, 50, 100 positions

Manual sampler

Dual furnace system Safety quick fit system Touch screen display

Standby mode

## **Physical Specifications**

Dimensions 81x50x37 cm

Weight 68 kg

Power supply 230V, 50/60Hz Adsorbed power 5A, 1100W

Gas requirements Helium (99.999% purity), 3-5 bar

Oxygen (99.999% purity), 3-5 bar

Air (oil free compressed air)



<sup>\*</sup> Accuracy and precision are related to samples nature and homogeneity .

## **Analytical and Technical Features**

## **Analytical Conditions**

Gas carrier Helium Leak test Automatic

Furnace temperature Left Furnace: max 1100°C RightFurnace: max 1100°C

Oven temperature Max 110°C

Oxygen volume need Automatically calculated by the oxygen doser

Flow rate Electronic Flow Rate
Gas separation 3 m GC Column
Detector High Sensitivity TCD

Software data analysis EAS Clarity

Calibration Linear, Quadratic, Cubic

Active calibration As needed

## **Sample**

Sample size Up to 1 g for soil samples

Sample type Liquid Solid

Capsule High purity tin and/or silver capsule

## **Accessories**

Installation kit
Sampling kits
Small samples; big samples; deluxe kit

Microbalance
Consumables
Proprietary NC Technologies S.r.I.

Technical assistance

By phone or email within 24 hours







Via Milano,15/A - 20060 Bussero (MI), Italy



Phone: +39 02 950 34 69



www.nctechnologies.it